## Autopsy Table ST 10/300 Fixed Mounted



Model ST 10/100 with additional tap beside the basin

The autopsy table is mounted on a center support column and is completely made of stainless steel (material no. 1.4301- ANSI 304). The work surface features a high, profiled surround and a seamlessly welded basin,  $400 \times 500 \times 200$  mm, is with a tube overflow and a sieve insert is installed at one end of the table. Additionally, the table top is equipped with a waterproof sound-proofing and features a negative inclination arranged diagonally to the drain. All power connections are built into the base and can be easily accessed by means of an inspection opening.

## The Autopsy Table ST 10/300 model includes the following features:

- Knee-operated mixing tap for cold and warm water at the center column
- 3 m shower hose with hand sprinkler
- · 2 splash proof electrical outlets
- · Electrical and water connections are preinstalled
- The safety requirements for drinking water as per DIN 1988 part 4 and EN DIN 1717 are to be provided by the customer

## **Technical Data:**

Dimensions W/D/H:  $2600 \times 850 \times 850 \text{ mm}$  Voltage: 230 Volt / 1 Ph / 50 Hz

 Carrying capacity:
 max. 300 kg

 Type:
 ST 10/300

 Article No:
 010.010.300

Alternatively available with a carrying capacity of 450 kg at an extra charge - Article No.: 010.010.858

KUGEL medical GmbH & Co. KG

Hermann-Köhl-Str. 2A

DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0 Fax +49/941/20 86 48-29 www.KUGEL-medical.de



