



CLASS I BIOSAFETY CABINET

CLASS I BIOSAFETY CABINET

Biolab Biological Safety Cabinet Class I series has an energy efficient design with low noise operation and corrosion free construction. Negative pressure ventilated cabinet operated with an open front and provides personnel protection. Air exiting the chamber is HEPA filtered to assure removal of contaminants.

Used in Biolab Biological Safety Cabinet Class I is widely used in microbiological work, research purposes in the field of genetics, cell biology, molecular biology, plant and animal cell culturing..

100 CLASS I BIOSAFETY CABINET



Stainless Steel 304 table for operation

HEPA Filter protects environment by filtering air before it is exhausted

UV lamp for sterilization

Soft touch control panel with LCD display

Stable air flow system providing personnel and environmental protection

Negative air pressure provides personnel protection by constant movement of air in working area

Low noise and high energy efficiency for operational cost savings

SPECIFICATIONS

Model	BCBS-101	BCBS-102	BCBS-103	BCBS-104
Type	Class I			
HEPA Filter	1			
HEPA Filter Efficiency	99.999% efficiency at 0.3 μ m			
Air Flow Velocity	> 0.3m/s	0.75 ~ 1.0 m/s	> 0.3m/s	0.75 ~ 1.0 m/s
Minimum Face Velocity	75 linear feet per minute (lfpm)			
LED Lamp	15Wx1	18Wx1	8Wx1	21Wx2
UV Lamp	15Wx1	20Wx1		
Work Surface Height	Benchtop Type			740 mm
Inner Dimension (WxDxH)	540x385x440 mm	680x540x500 mm	800x635x600 mm	968x695x630 mm
Overall Dimension (WxDxH)	550x395x730 mm	700x550x900 mm	900x713x1250 mm	1100x695x1924 mm
Package Size	710x550x950 mm	870x720x1170 mm	1060x870x1440 mm	1270x880x2230 mm
Chamber Material	Poly Methyl Methacrylate (PMMA)		Cold-rolled steel coated with anti-bacteria powder	
Noise Level	\leq 60 dB	\leq 55 dB	\leq 60 dB	\leq 67 dB
Weight	40 kg	60 kg	100 kg	150 kg
Power	150 W	\leq 160 W	\leq 400 W	\leq 350 W
Power Supply	AC220V \pm 10%,50/60Hz;110V \pm 10%,60Hz			
Caster	-	-	-	Universal Wheel

ACCESSORIES

Accessory Code	Name	For Models
BX-103986	Fluorescent Lamp	All
BX-103987	UV Lampx2	All

BX-103988	Base Stand	BCBS-104
-----------	------------	----------

OPTIONAL ACCESSORIES

Accessory Code	Name	External Dimension	Package Size (mm)
BX-102683	Infrared Sterilizer	150x95x210 mm	
BX-102684	Airflow Tester		
BX-102685	Formalin Fumigation Sterilizer		300x200x160 mm
BX-102686	Ammonium Hydrogen Carbonate Neutralizer		300x200x160 mm
BX-102687	Armrest		
BX-102688	Laboratory Chair		
BX-102689	Dust Particle Counter		
BX-102690	Digital Sound Level Meter		
BX-102691	Illumination Meter		
BX-102692	Air Flow Anemometer		
BX-102693	Protective Garment		
BX-102694	Protective Gloves		



BCBS-101, BCBS-102



BCBS-103



BCBS-104



Biolab Scientific Limited

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: contact@biolabscientific.com | Phone: +1 416 639 2438 | Website: www.biolabscientific.com