

Down-Draft Grossing Table UCS-E-1500-D

Electrical height adjustable – with Illumination



Table with accessories (not included in delivery)

The down-draft grossing table is equipped with air down draft system. This specialized work table guarantees a laboratory environment free of pollutants while slicing and preparing histological slide preparations and is even well suited for small dissections. Airflow travels over the entire working area and is vacuumed out downwards, which provides optimal removal of unsanitary vapours. MAC (maximum allowable concentration) values for formalin are below target throughout the entire vacuum-effected area. The cutting table is to be connected to an exhaust system provided by the customer. The UCS-E-1500 with illumination includes the following features:

Workstations with a large surface area

- Entire construction is made of stainless steel (Material no. 1.4301 – ANSI 304)
- Stable, self-supporting undercarriage with adjustable legs
- Low overall height of the upper structure increases leg room
- Height adjustable from 750 to 1000 mm
- Height adjustable with UP/Down push bottom
- Work surface includes a seamlessly welded basin measuring 400 x 400 x 200 mm
- Basin includes integrated air vacuum
- Very sturdy, perforated and insertable working plates
- Integrated sprinkler system with adjustable valve for bath sanitation located underneath the perforated working plates
- Formalin sink with sieve insert
- Mixing tap with detachable cold and hot water shower
- Exhaust duct
- Build in Illumination with “Daylight” Quality approx. 5000 °K for direct photograph without additional Flashlights
- 2x 230 Volt Electric socket outlet

Technical data see back page

KUGEL Medizintechnik Vertriebs GmbH
Hermann-Köhl-Str. 2A
DE – 93049 Regensburg / Germany

Phone +49/941/208648-0
Fax +49/941/208648-29
www.KUGEL-medical.de



KUGEL
medical
■■■■

Down-Draft Grossing Table UCS-E-1500-D

Electrical height adjustable – with Illumination

Technische Details:

Dimensions Table W/D/H:	1500 x 1000 x 750 / 1000 mm
Dimensions Illumination:	1500 x 300 x 1100 mm
Exhaust Air capacity:	Min. 1000 m³/h
Material:	Geschliffenes Edelstahl (1.4301 – ANSI 304)
Perforated working area:	ca. 960 x 920 mm
Cold / Warm water connection:	DN 15
Waste water connection:	DN 50
Exhaust duct diameter:	DN 250
Type:	UCS-E-1500-D-BEA
Article-No.:	060.020.100

Accessories:



Waste bin on the UCS
Art.-No. 060.012.056



Magnifier with illumination
Art.-No. 060.015.050



Kameraholder on the railingsystem KA-05
Art.-No. 060.015.081



Dictaphone holder Art.-No. 060.020.041
Document holder Art.-No. 060.012.031
Formalin dispensing system Art.-No. 060.012.024



Cutting board "CORIAN"
with ruler Art.-No. 060.019.030

KUGEL Medizintechnik Vertriebs GmbH
Hermann-Köhl-Str. 2A DE – 93049 Regensburg/Germany

www.KUGEL-medical.de
Phone +49/9 41/20 86 48-0
Fax +49/9 41/20 86 48-29

UNI BRAND INTERNATIONAL
Flat No 104, Mohammadia Plaza,
College Road, Rawalpindi Pakistan.
Tel : +92-51-5534794, 5541767, 5557602
Fax : +92-51-5555392
Cell : +92-345-5577591
Lahore (Cell): +92-345-8568214
Rawalpindi Lahore Karachi



KUGEL
medical