

Fine-Tek Plastic Processing/Embedding Cassettes

Biolab USA



Fine-Tek Plastic Processing/Embedding Cassettes

FINE-Tek PLASTIC today can manufacture and supply many models of cassettes for processing and embedding histological tissue specimens, most of which can be adapt to different automated cassettes labeling system.

FINE-Tek PLASTIC EM 105 Biopsy Processing/Embedding Cassettes

Special hinge design locks the cassette and lid together firmly, and realizes more frequent opening & closing of the lid without any danger of specimen loss

- Made of P.O.M material, 100% resistant to chemical reaction of histological solvents
- With square holes in dimensions of 0.9mm×0.9mm
- With two large labeling areas on the cassette sides and one slanted front writing surface of 45° angle, adaptable to most of the cassette labelers
- Detachable lid easy to remove and re-insert
- pre-loaded by sleeves or tapes, considerably optimizing the user convenience

COLOR	UNASSEMBLED PACKED BY DISPENSER BOX	ASSEMBLED PACKED BY DISPENSER BOX	PRE-LOADED IN SLEEVES	TAPPED
White	0103-1100-16	0103-1102-16	0103-1103-16	0103-1132-16
Yellow	0103-1100-17	0103-1102-17	0103-1103-17	0103-1132-17
Green	0103-1100-08	0103-1102-08	0103-1103-08	0103-1132-08
Pink	0103-1100-12	0103-1102-12	0103-1103-12	0103-1132-12
Blue	0103-1100-02	0103-1102-02	0103-1103-02	0103-1132-02
Tan	0103-1100-15	0103-1102-15	0103-1103-15	0103-1132-15
Purple	0103-1100-09	0103-1102-09	0103-1103-09	0103-1132-09
Peach	0103-1100-11	0103-1102-11	0103-1103-11	0103-1132-11
Orange	0103-1100-10	0103-1102-10	0103-1103-10	0103-1132-10
Aqua	0103-1100-00	0103-1102-00	0103-1103-00	0103-1132-00

Note: We have all types of Plastic Processing/Embedding Cassettes.