



MICROSCOPE COVER GLASS

SIZE 22X22 & 24X50



Non-Sticking Microscope Cover Glass

NEW

-ASM surface treatment

Borosilicate glass 3.3 material, with special surface treatment, especially suitable for automatic sealing machine

- High cleanliness, free cleaning, the surface is specially treated, and it is not easy to stick under normal storage conditions;

- Appropriate flexibility, withstand a certain degree of bending without breaking;

- The surface characteristics make the sealant evenly distributed to prevent air bubbles from entering, and the specimens can be stored for a long time.

Borosilicate 3.3 Glass

SPECIFICATION	PACKING
24×40mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×50mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×60mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×40mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×50mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×60mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton

Super White Glass

SPECIFICATION	PACKING
24×40mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×50mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×60mm, #1.0 thickness 0.13-0.16mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×40mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×50mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton
24×60mm, #1.5 thickness 0.16-0.19mm	200pcs/box, 1,000pcs/innerbox, 20,000pcs/carton

