



LASER CASSETTE PRINTER Fine Print FP – 600



Character:

- Simple design, easy use;
- High printing speed, 2.5-4 seconds per cassette;
- Non-contact UV laser technology;
- Excellent ability to export characters, numbers, letters, 1D & 2D codes, graphics;
- Fully compatible with LIS/HIS system;
- Support color setting print;
- 10 inch LED touch screen via built-in software, also support to connect external PC;
- Built in bar-code scanner, and support to connect external scanners;
- Built in air purification system ;
- Compatible with various models of embedding cassettes on the markets.

Technical Data:

Printing technology	UV Laser marking
Print speed	2.5-4 seconds/cassette
Print resolution	2500 dpi
Cassette type	Cassette with lid assembled(two-piece); Cassette without lid
Loading hopper	6
Max Capacity Per Loading hopper	100 cassettes with lid assembled or without lid
Unload station max capacity	119 cassettes
Compatible system	Open to connect HIS/LIS systems
Dimension LxWxH	486mmx313mmx497mm
Weight	37kg
Power requirement	100-240 V AC, 50/60 Hz, 300W
Other configuration	1. 10 inch touch screen 2. Built-in bar-code scanner; 3. Built in air purification system.

