



TISSUE EMBEDDING CENTER

Micro Tek 6LA



Features:

- The newest high technology and Button Panel Control,new type of heating material.
- The machine is made up of embedding machine; freezing part .the design can be flexibly selected between embedding machine and freezing part.
- Automated processor control systems make it possible for presetting the start time and temperature of any day.
- special nano technical make it heat up quickly and equably.

- The LCD display can show you the temperatures of the paraffin bath, thermostatic bath and paraffin nozzle. So the user can operator it easily.
- inductive head lamp and metal switch, user can touch the switch by tweezers or knives to control the paraffin.
- Automatic memory the presetting temperatures and start to work based on the time.
- large capacity crock can contain enough paraffin for use at one time.
- Multilevel security protection system.
- There are 3 heating hole for forceps in two side of the paraffin nozzle.

Specification:

- Capacity of the wax bath : 6000 ml.
- Melting paraffin prefabricate: 0 ~99°C
- Freezing range temperature: 0 ~ -20°C
- Complete machine consuming more electricity: Less than 1KVA
- Control mode: Auto, foot
- Environment temperature: 0°C~ 40°C
- Time control: Opening and closing time of the machine is set in 24 hours optionally.
- nominal voltage: AC220V±10%.
- rated frequency: 50±1Hz
- Dimension: 1180×660×470mm
- weight: 87kg