

Micro Tek BL – 100 Pathology Biopsy Table Grossing Work Station



Dimensions		150 ~ 200cm length width × 190 × 80cm high
Air suction		≥ 40 cubic meters / per minute
Noise		≤50db
Features		
1	All stainless steel pull pattern	Constitutes high-quality steel, corrosion resistant
2	Air and suction systems	High-power suction devices, low noise. Quickly to remove harmful gases.
3	Built-in flush device table	Can quickly remove water
4	Fast pressure liquid sewage discharge	To completely remove water
5	With pipe mill	Organization of debris to avoid blocking pipeline
6	Configuration disinfection lamps, lighting	

7	Magnetic tool holder adsorption	Easy to use
8	With warm water	With warm water, easy winter drawn
9	Rinse screen with a small sample	Easy cleaning of small specimens to avoid the loss of small specimens
10	With sink	
11	Optional large specimen digital imaging system	Improve efficiency, reduce errors and data permanently archived for easy access. Gross specimen images simultaneously for reference, diagnosis is based on more fully