

Post Mortem Instruments 11-975 / Set 2



Post Mortem Instrument Set 2 consisting of:

- 1 pc. Autopsy knife with metal handle, cutting length 195 mm
 - 1 pc. Dissecting knife, cutting length 100 mm
 - 1 pc. Autopsy knife with wooden handle, cutting length 110 mm
 - 1 pc. Cartilage knife with wooden handle, Cutting length 70 mm
 - 2 pc. Incision scissors, 145 mm
 - 1 pc. Bowel scissors, 210 mm
 - 1 pc. Dissecting forceps, standard, length 160 mm
 - 1 pc. Dissecting forceps, 2 x 3 teeth, length 145 mm
 - 1 pc. Dissecting forceps 1 x 2 teeth, length 250 mm
 - 1 pc. Directors 160 mm
 - 2 pc. Probes 250 mm / 1,5 mm
 - 1 pc. Post mortem needles, 115 mm (12 pc.)
 - 1 pc. Folding ruler 300 mm
 - 1 pc. Tape measure in plastic capsule, 1500 mm
 - 1 pc. Bow Saw, 420 mm incl. Spare saw blade
 - 1 pc Bone shears 230 mm
 - 1 pc. Bone holding forceps, Langenbeck 215 mm
 - 1 pc. Chisel, 245 mm
 - 1 pc. Metal mallet, Ombredanne, 240 mm, 550 gram
 - 1 pc. Skull breaker 140 mm
 - 1 pc. Microscopic needle holder with plastic handle, 240 mm
 - 1 pc. Dermatograph, red
- Article-No.: 020.010.196

Alternatively:

Post Mortem Instrument Set 2 in aluminium case
Article-No.: 020.010.405

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■