

SIGNATURE SLIDE PRINTER

USA, Canada, Latin America:
1-800-797-2772 or +1-763-475-6676
sales@primerahealthcare.com

Europe, Middle East, Africa:
+49-611-927770
sales@primera.eu

Asia Pacific:
+61 3 9798 8720
sales@primera-ap.com

Print directly onto slides for more efficient and more accurate specimen identification.

Primera's new Signature Slide Printer can significantly increase the efficiency of your lab while helping to reduce the risk of misidentification of specimens. It prints directly onto slides, eliminating handwriting or expensive, difficult to apply xylene-resistant labels. With impressive 300 dpi print resolution, you can print text, graphics and logos along with high-resolution linear and 2D bar codes on every slide you process.



Key Features and Benefits:

- ✓ On-demand, full-color printing – print only the number of slides you need to help eliminate errors
- ✓ Reduce inventory and purchasing activity – your inventory of multiple colors of slides is eliminated since you can now print colors on white slides
- ✓ Up to 9 slides per minute
- ✓ Output capacity: 15 slides
- ✓ Prints directly onto slides – eliminates handwriting that is hard to read and labels that are hard to apply and can potentially fall off
- ✓ Small footprint – fits perfectly next to the microtome
- ✓ Chemical-resistant ink – ensures reliable identification of slides
- ✓ USB interface – ability to integrate with Laboratory Information Systems

Easy to Use

Signature is designed for efficient, hands-free operation. Slides are stored in easy-to-load cartridges – away from dust and other potential contaminants. It takes just a few seconds to change out the cartridge if you need standard slides for one study and positively-charged slides for others. The level of slides remaining in a cartridge is easily viewed through the cartridge's transparent blue LED backlit case.

Optional PTLab™ Software

PTLab Software was designed by Primera specifically for use with the Signature Slide Printer. It can be installed on any Windows Vista 7/8/10 PC. With PTLab, you'll be up and running and printing full-color slides in just minutes. Printing only the slides you need, when and where you need them, helps avoid potentially serious patient safety errors.

