

Tissue-Tek Prisma® Plus & Film®

Automated Slide Stainer & Coverslipper



SMART Automation
Powered by Sakura

Specifications

Capacity	660 slides
Dimensions	
W x D x H	197 x 71 x 69 cm
Weight	225 kg
Capacity	
Loading	Up to 60 slides
Unloading	Up to 240 slides

Specifications Prisma® Plus Slide stainer

Capacity	660 slides
Dimensions	
W x D x H	125 x 71 x 63 cm
Weight	150 kg
Number of reservoirs	
Start	1 to 3 stations
Reagent	31 to 51 reservoirs
Drying	0 to 2 stations
Wash	0 to 4 stations (8 programmable steps)
Heating (optional)	2 stations
Capacity	
Reagent volume	160, 260 and 680 mL
Staining protocols	11 runs continuously loading, and simultaneously running
Baskets	10- or 20-slide basket
Loading	Up to 60 slides per run
Temperature	
Drying station	30 ~ 65°C
Heating station	30 ~ 70°C

Specifications Tissue-Tek Film® Coverslipper

Capacity	1080 slides
Dimensions	
W x D x H	72 x 58 x 69 cm
Weight	75 kg
Capacity	
Loading	Up to 60 slides
Unloading	Up to 240 slides
Track & Trace	Slide Barcode Reading

Tissue-Tek Prisma® Plus & Tissue-Tek Film®

The Tissue-Tek Prisma® Plus & Tissue-Tek Film® integrated system is the gold standard in staining and film coverslipping. This reliable and flexible workhorse is the ultimate timesaver for your laboratory.

This guarantees true walk away capability and continuous confidence in performance and quality due to its unsurpassed system up-time.

The flexibility of standard and special staining in a single system, with the use of proven film coverslipping technology makes this a cost efficient method.

Ultimate time saving

Virtually every histopathology lab is pushed on turnaround times and faces increasing workload, often with reduced staff and with equal or less financial resources. The integrated Prisma Plus and Film is the ultimate time saver, freeing your staff for other important tasks.

Supreme flexibility

Customers have greatly appreciated the convenience delivered by the integrated Prisma Film stainer coverslipper. This is due to the flexibility of having a single platform to easily perform a H&E and special stains reducing hands on time, maintenance time and service costs.

Superior reliability

The Prisma Plus and Film coverslipper have been designed, engineered and built to provide the laboratory with an integrated stainer coverslipper having a pedigree of dependability and reliability. This ensures a smooth workflow through continuous loading and high quality staining slide after slide.

Tissue-Tek Prisma Plus & Tissue-Tek Film offer you:

- Ultimate time saving
- Supreme flexibility
- Superior reliability



Tissue-Tek® iSupport™ Powered by Sakura

Item code Description

8022	License Tissue-Tek® iSupport™
8023	License Tissue-Tek® iSupport™ (extra user)

Item code Description

6172	Tissue-Tek Prisma® Plus Standard Configuration; 230 V, 50/60 Hz
6173	Tissue-Tek Prisma® Plus Special Staining Configuration; 230 V, 50/60 Hz
4768	Tissue-Tek® 20-Slide Basket
6135	Tissue-Tek Prisma® 20-Slide Basket Load/Unload Adapter
6136	Tissue-Tek Prisma® 20-Slide Basket Adapter
6137	Tissue-Tek Prisma® 10-Slide Basket
6138	Tissue-Tek Prisma® 10-Slide Basket Adapter
6139	Tissue-Tek Prisma® 10-Slide Basket Load/Unload Station Adapter
6140	Tissue-Tek Prisma® Special Stain Solution Reservoir (160 mL)
6141	Tissue-Tek Prisma® Heated Solution Reservoir
6144	Tissue-Tek Prisma® Lid for Special Stain Solution Reservoir
6145	Tissue-Tek Prisma® Small Solution Reservoir (260 mL)
6146	Tissue-Tek Prisma® Lid for Small Solution Reservoir
6147	Tissue-Tek Prisma® Standard Solution Reservoir (680 mL)
6148	Tissue-Tek Prisma® Standard Solution Reservoir with Handle (680 mL)
6149	Tissue-Tek Prisma® Wash Reservoir
6151	Tissue-Tek Prisma® Lid for Standard Solution Reservoir
6152	Tissue-Tek Prisma® 3-Position Reservoir Tray
6153	Tissue-Tek Prisma® Large Lid for 3-Position Reservoir Tray
6154	Tissue-Tek Prisma® 4-Position Reservoir Tray
6155	Tissue-Tek Prisma® Large Lid for 4-Position Reservoir Tray
6156	Tissue-Tek Prisma® Special Stain Reservoir Tray
6157	Tissue-Tek Prisma® Standard Solution Reservoir Set
6158	Tissue-Tek Prisma® Small Solution Reservoir Set
6159	Tissue-Tek Prisma® Special Stain Solution Reservoir Set



Staining & Coverslipping

Item code Description

6160	Tissue-Tek Prisma® Activated Carbon Filter; 2/pcs
6162	Tissue-Tek Prisma® Basket Adapter Label
6163	Tissue-Tek Prisma® Control Panel Protection Sheet
6164	Tissue-Tek Prisma® Compact Flash Card (512 Mb)
6165	Tissue-Tek Prisma® Duct Connection Adapter (38 mm)
6166	Tissue-Tek Prisma® Duct Connection Adapter (75 mm)
HIS40614	Tissue-Tek® Insert for slide basket (for 5 large slides maximum 60 mm wide)

SA-6172-UT	User Training Tissue-Tek Prisma® Plus Standard
SA-6173-UT	User Training Tissue-Tek Prisma® Plus Special
SC-6172-G	Service Agreement Tissue-Tek Prisma® Plus Standard Gold
SC-6172-P	Service Agreement Tissue-Tek Prisma® Plus Standard Platinum
SC-6172-S	Service Agreement Tissue-Tek Prisma® Plus Standard Silver
SC-6173-G	Service Agreement Tissue-Tek Prisma® Plus Special Gold
SC-6173-P	Service Agreement Tissue-Tek Prisma® Plus Special Platinum
SC-6173-S	Service Agreement Tissue-Tek Prisma® Plus Special Silver

Item code Description

4742	Tissue-Tek Film® Coverslipper 230 V, 50/60 Hz
6508	Tissue-Tek® Barcode system for Coverslippers
6134	Tissue-Tek Prisma® Link System Kit for Tissue-Tek Film® Coverslipper
4770	Tissue-Tek® Coverslipping Film 70 m
4772	Tissue-Tek® Coverslipping Film Cutter Blades; 5/cs
6160	Tissue-Tek Prisma® Activated Carbon Filter; 2/cs
4745	Tissue-Tek® Xylene Bottle (500 mL)
4746	Tissue-Tek® Cap for Xylene Bottle
4748	Tissue-Tek® Waste Bottle (200 mL)
4749	Tissue-Tek® Cap for Waste Bottle
6506	Tissue-Tek® Exhaust hose assembly 38 mm Ø
6507	Tissue-Tek® Exhaust hose assembly 75 mm Ø

SA-4742-UT	User Training Tissue-Tek Film® Coverslipper
SC-4742-G	Service Agreements Tissue-Tek Film® Coverslipper Gold
SC-4742-P	Service Agreements Tissue-Tek Film® Coverslipper Platinum
SC-4742-S	Service Agreements Tissue-Tek Film® Coverslipper Silver