

Tissue-Tek Prisma® Plus & Tissue-Tek® Glas™ g2 Automated Slide Stainer & Coverslipper



SMART Automation
Powered by Sakura

Specifications

Throughput	360 slides per hour
Dimensions	
W x D x H	200 x 71 x 75 cm
Weight	260 kg
Capacity	
Loading	Up to 60 slides
Unloading	Up to 240 slides

Specifications Prisma® Plus Slide stainer

Capacity	660 slides
Dimensions	
W x D x H	125 x 71 x 63 cm
Weight	150 kg
Number of reservoirs	
Start	1 to 3 stations
Reagent	31 to 51 reservoirs
Drying	0 to 2 stations
Wash	0 to 4 stations (8 programmable steps)
Heating (optional)	2 stations
Capacity	
Reagent volume	160, 260 and 680 mL
Staining protocols	11 runs continuously loading, and simultaneously running
Baskets	10- or 20-slide basket
Loading	Up to 60 slides per run
Temperature	
Drying station	30 ~ 65°C
Heating station	30 ~ 70°C

Specifications Tissue-Tek® Glas™ g2 Coverslipper

Throughput	360 slides per hour
Dimensions	
W x D x H	75 x 62 x 75 cm
Weight	110 kg
Capacity	
Loading	Up to 60 slides
Unloading	Up to 240 slides

Tissue-Tek Prisma® Plus & Tissue-Tek® Glas™ g2

Tissue-Tek Prisma® Plus & Tissue-Tek® Glas™ g2 provides high reliability, user convenience and flexibility in staining and glass coverslipping. One system for routine and special staining, with Integrated Track & Trace with barcode reading.

Optimised consistent high quality

Prisma Plus & Tissue-Tek Glas g2 enables optimised and consistent high quality routine and special staining. Slide-to-slide consistency results in crisp images with clear morphological and nuclear detail. Quality is further enhanced with optimised protocols and specialised Sakura application support. Integrated track & trace capability through slide barcode reading enhances laboratory quality assurance. Determine the status of slides real-time and stay in control.

Increased efficiency and reliability

Prisma & Glas g2 provides unmatched staining and glass coverslipping turn around time improvement with 360 slides per hour. The system enables an efficient workflow with continuous loading and unloading capabilities. The Tissue-Tek Prisma Plus is the most flexible multiple stainer on the market, with a wide range of configurations for routine and special staining. The availability of three different volume reagent reservoirs optimises reagent usage and decreases operating costs. Mechatronics™ mechanism prevents any damage to the cover glass or tissue and guarantees bubble free slides.

Enhance user convenience with Tissue-Tek® Glas™ g2 Coverslips & Easy Holders. These pre-filled coverslip holders enable quick and easy refilling of glass coverslips.

Tissue-Tek Prisma & Glas g2 offer you:

- Improved quality, efficiency and user friendliness
- Rapid glass coverslipping with 360 slides per hour
- One system for routine and special staining
- Integrated Track & Trace and slide identification with barcode reading



Item code Description

6172	Tissue-Tek Prisma® <i>Plus</i> Standard Configuration; 230 V, 50/60 Hz
6173	Tissue-Tek Prisma® <i>Plus</i> Special Staining Configuration; 230 V, 50/60 Hz
4768	Tissue-Tek® 20-Slide Basket
6135	Tissue-Tek Prisma® 20-Slide Basket Load/Unload Adapter
6136	Tissue-Tek Prisma® 20-Slide Basket Adapter
6137	Tissue-Tek Prisma® 10-Slide Basket
6138	Tissue-Tek Prisma® 10-Slide Basket Adapter
6139	Tissue-Tek Prisma® 10-Slide Basket Load/Unload Station Adapter
6140	Tissue-Tek Prisma® Special Stain Solution Reservoir (160 mL)
6141	Tissue-Tek Prisma® Heated Solution Reservoir
6144	Tissue-Tek Prisma® Lid for Special Stain Solution Reservoir
6145	Tissue-Tek Prisma® Small Solution Reservoir (260 mL)
6146	Tissue-Tek Prisma® Lid for Small Solution Reservoir
6147	Tissue-Tek Prisma® Standard Solution Reservoir (680 mL)
6148	Tissue-Tek Prisma® Standard Solution Reservoir with Handle (680 mL)
6149	Tissue-Tek Prisma® Wash Reservoir
6151	Tissue-Tek Prisma® Lid for Standard Solution Reservoir
6152	Tissue-Tek Prisma® 3-Position Reservoir Tray
6153	Tissue-Tek Prisma® Large Lid for 3-Position Reservoir Tray
6154	Tissue-Tek Prisma® 4-Position Reservoir Tray
6155	Tissue-Tek Prisma® Large Lid for 4-Position Reservoir Tray
6156	Tissue-Tek Prisma® Special Stain Reservoir Tray
6157	Tissue-Tek Prisma® Standard Solution Reservoir Set
6158	Small Solution Reservoir Set
6159	Special Stain Solution Reservoir Set
6160	Tissue-Tek Prisma® Activated Carbon Filter; 2/cs
6162	Tissue-Tek Prisma® Basket Adapter Label
6163	Tissue-Tek Prisma® Control Panel Protection Sheet
6164	Tissue-Tek Prisma® Compact Flash Card (512 Mb)
HIS40614	Tissue-Tek Prisma® Insert for slide basket (for 5 large slides maximum 60 mm wide)



Item code Description

SA-6172-UT	User Training Tissue-Tek Prisma® Plus Standard
SA-6173-UT	User Training Tissue-Tek Prisma® Plus Special
SC-6172-G	Service Agreement Tissue-Tek Prisma® Plus Standard Gold
SC-6172-P	Service Agreement Tissue-Tek Prisma® Plus Standard Platinum
SC-6172-S	Service Agreement Tissue-Tek Prisma® Plus Standard Silver
SC-6173-G	Service Agreement Tissue-Tek Prisma® Plus Special Gold
SC-6173-P	Service Agreement Tissue-Tek Prisma® Plus Special Platinum
SC-6173-S	Service Agreement Tissue-Tek Prisma® Plus Special Silver

Item code Description

6502	Tissue-Tek® Glas™ g2 Coverslipper 230V, 50/60Hz
6508	Tissue-Tek® Barcode System for Coverslippers
6168	Tissue-Tek Prisma® Link System Kit for Tissue-Tek® Glas™ g2 Coverslipper
6415	Tissue-Tek® Glas™ g2 Coverslips & Easy Holder 24x40 mm; 5 x 200/cs
6416	Tissue-Tek® Glas™ g2 Coverslips & Easy Holder 24x50 mm; 5 x 200/cs
6417	Tissue-Tek® Glas™ g2 Coverslips & Easy Holder 24x60 mm; 5 x 200/cs
6423	Tissue-Tek® Glas™ g2 Nozzle reservoir
6430	Tissue-Tek Prisma® Waste Container
6433	Tissue-Tek Prisma® bottle
6504	Tissue-Tek® Glas™ g2 Slide Rack
6505	Tissue-Tek® Glas™ g2 Waste bottle
6160	Tissue-Tek Prisma® Activated Carbon Filter; 2/cs
6506	Tissue-Tek® Exhaust hose assembly 38 mm Ø
6507	Tissue-Tek® Exhaust hose assembly 75 mm Ø
SA-4742/	User Training Tissue-Tek® Glas™ g2
6502-UT	Coverslipper
SC-6502-G	Service Agreements Tissue-Tek® Glas™ g2 Coverslipper Gold
SC-6502-P	Service Agreements Tissue-Tek® Glas™ g2 Coverslipper Platinum
SC-6502-S	Service Agreements Tissue-Tek® Glas™ g2 Coverslipper Silver

Staining & Coverslipping

Tissue-Tek® iSupport™ Powered by Sakura

Item code Description

8022	License Tissue-Tek® iSupport™
8023	License Tissue-Tek® iSupport™ (extra user)