



Vacumm Blood Collection Tube

Finelabware

Vacumm Blood Collection Tube

A Vacutainer blood collection tube is a sterile glass or plastic test tube with a colored rubber stopper creating a vacuum seal inside of the tube, facilitating the drawing of a predetermined volume of liquid. Analysis of blood samples has been one of the most important diagnostic tools available for the detection and treatment of diseases. Although blood samples in minute quantities can be taken from finger tips.

Vacumm Blood collection tube Plain Tube (No additive)

A	B	C	D	E
Material: glass	Material: glass	Material: PET	Material: PET	Material: PET
Capacity: 1-5ml	Capacity: 6-7ml	Capacity: 1-5ml	Capacity: 6-7ml	Capacity: 8-10ml
Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 16×100mm
Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray
18trays/ctn	18trays/ctn	18trays/ctn	18trays/ctn	12trays/ctn

Vacumm Blood collection tube Clot activator tube

A	B	C	D	E
Material: glass	Material: glass	Material: PET	Material : PET	Material: PET
Capacity: 1-5ml	Capacity: 6-7ml	Capacity: 1-5ml	Capacity : 6-7ml	Capacity: 8-10ml
Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 16×100mm
Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray
18trays/ctn	18trays/ctn	18trays/ctn	18trays/ctn	12trays/ctn

Vacumm Blood collection tube EDTA K2/K3

A	B	C	D	E
Material: glass	Material: glass	Material: PET	Material: PET	Material: PET
Capacity: 1-5ml	Capacity: 6-7ml	Capacity: 1-5ml	Capacity: 6-7ml	Capacity: 8-10ml
Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 13×75mm	Spec.: 13×100mm	Spec.: 16×100mm
Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray	Package: 100pcs/tray
18trays/ctn	18trays/ctn	18trays/ctn	18trays/ctn	12trays/ctn

Vacuum Blood collection tube ESR tube (sodium citrate 1:4)

A	B	C	D	E
Material:glass Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:glass Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:8-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

Vacuum Blood collection tube gel+Clot activator

A	B	C	D	E
Material:glass Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:glass Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:1-5ml Spec.:13×75mm Package:100pcs/ tray 18trays/ctn	Material:PET Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:8-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

Non vacuum blood collection tube Plain tube

A	B	C
Material:PP Capacity:1-4ml Spec.:13×75mm Package:100pcs/tray 12trays/ctn	Material:PP Capacity:5-7ml Spec.:13×100mm Package:100pcs/tray 12trays/ctn	Material:PET Capacity:7-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

Non vacuum blood collection tube Micro tub

A
Material:PP Capacity:0.5ml Spec.:8×45mm Package:100pcs/tray 30trays/ctn

Vacuum Blood collection tube Heparin sodium (Heparin lithium)

A	B	C	D	E
Material:glass Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:glass Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:8-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

Vacuum Blood collection tube PT tube (sodium citrate 1:9)

A	B	C	D	E
Material:glass Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:glass Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:8-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

Vacuum Blood collection tube Glucose (sodium fluoride/potassium oxalate)

A	B	C	D	E
Material:glass Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:glass Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:1-5ml Spec.:13×75mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:6-7ml Spec.:13×100mm Package:100pcs/tray 18trays/ctn	Material:PET Capacity:8-10ml Spec.:16×100mm Package:100pcs/tray 12trays/ctn

