

-25°C FREEZER



-25°C FREEZER

A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency. Safety features and durable interiors provide you with the most dependable performance.

100 SERIES -25°C UPRIGHT FREEZER



- Microprocessor temperature controller
- LCD display
- ABS liner interior and painted steel board exterior
- Door with safety lock and rotatable sure-grip door handle
- Polycarbonate insulated inner doors
- Insulation thickness: 90 mm
- Refrigerant: CFC free R507
- Password protected configuration page
- Danfoss compressor and EBM fan motor
- 7 drawers made of plastic inside
- 4 casters for easy handling
- Safety functions: Audible/visual alarm for high or low temperature and system failure

SPECIFICATIONS

Model	BFUR-25-101
Capacity	270 L
Temperature Range	-10°C to -25°C
Defrosting	Manual
Internal Dimension	440x440x1490 mm
Overall Dimension	580x660x1750mm
Package Size	655x695x1820 mm
Drawers/Shelves	7 drawers
Weight (Net/Gross)	65/70 kg
Input Power	370 W
Power Supply	220V,50/ 60Hz, 110V,50/60Hz

200 SERIES -25°C FREEZER CHEST



Stainless Steel exterior and interior made of aluminum
 Digital temperature display
 Rounded corners door design with lock
 Polycarbonate insulated inner doors
 Insulation thickness: 40mm
 Quick freezing function, i.e., 2-3 hours for -25°C
 High quality compressor: Low energy consumption, low noise level and low heat dissipation
 Refrigerant: CFC free R134a
 Supplied with Steel Wire Basket (1 unit) to store samples
 4 casters for easy handling
 Safety functions: Audible/visual alarm for high or low temperature alarm, system failure alarm

-25°C FREEZER

SPECIFICATIONS

Model	BFRZ-25-201	BFRZ-25-202	BFRZ-25-203
Capacity	110 L	166 L	196 L
Temperature Range	-10°C to -25°C		
Defrosting	Manual		
Overall Dimension	550x585x855 mm	610x875x840 mm	610x1075x768 mm
Package Size	625x635x900 mm	615x925x900 mm	615x1130x925 mm
Weight (Net/Gross)	35/45 kg	38/48 kg	50/60 kg
Input Power	145 W	160 W	180 W
Power Supply	220V, 50Hz		

Model	BFRZ-25-204	BFRZ-25-205	BFRZ-25-206
Capacity	226 L	358 L	508 L
Temperature Range	-10°C to -25°C		
Defrosting	Manual		
Overall Dimension	610x1072x840mm	670x1295x908 mm	840x1590x862mm
Package Size	754x1130x990 mm	855x1645x1070 mm	
Weight (Net/Gross)	55 kg	80/90 kg	100/110 kg
Input Power	207 W	320 W	375 W
Power Supply	220V, 50Hz		



Biolab Scientific Limited

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: contact@biolabscientific.com | Phone: +1 707 533 1445 | Website: www.biolabscientific.com