

-40°C FREEZER



-40°C FREEZER

Equipped with precise temperature controls, rugged interiors and safety features, our freezers provide you with the most dependable performance with low maintenance. A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency.

100 SERIES -40°C UPRIGHT FREEZER



- Microprocessor temperature controller
- LCD display
- Danfoss compressor and EBM fan motor
- Exterior made of painted steel board and interior of ABS liner
- Door with safety lock and rotatable sure-grip door handle
- Refrigerant: CFC free R507
- 4-7 drawers made of plastic inside
- 4 casters for easy handling
- Safety functions: Audible/visual alarm for high or low temperature and system failure

SPECIFICATIONS

Model	BFUR-40-101	BFUR-40-103	BFUR-40-104
Capacity	90 L	270 L	362 L
Temperature Range	-10°C to -40°C		
Defrosting	Manual		
Internal Dimension	438x388x539 mm	440x400x1490 mm	658x426x1290 mm
Overall Dimension	623x544x842 mm	570x595x1760 mm	930x731x1627 mm
Drawers/Shelves	4 drawers	6 drawers	4 drawers
Weight (Net/Gross)	44/48 kg	65/70 kg	102/110 kg
Input Power	280 W	300 W	400 W
Power Supply	220V, 50Hz	220V,50/ 60Hz, 110V,50/60Hz	220V,50/ 60Hz, 110V,50/60Hz
Package Size	-	650x750x1825 mm	820x1040x1850 mm



BFUR-40-101



BFUR-40-103



BFUR-40-10

200 SERIES -40°C FREEZER CHEST



Microprocessor temperature control
Digital temperature display
Stainless steel interior and painted steel board exterior
Door with safety lock
Polycarbonate insulated inner doors
Insulation thickness: 90 mm
Refrigerant: CFC free R507
High quality compressor and EBM fan motor: Low energy consumption, low noise level and low heat dissipation
Supplied with Steel Wire Basket (1 unit) to store samples
4 casters for easy handling
Safety functions: Audible/visual alarm for high or low temperature alarm, system failure alarm

-40°C FREEZER

SPECIFICATIONS

Model	BFRZ-40-201	BFRZ-40-202	BFRZ-40-203
Capacity	110 L	251 L	351 L
Temperature Range	-10°C to -40°C		
Defrosting	Manual		
Internal Dimension	410x410x654 mm	430x1050x556 mm	-
Overall Dimension	549x549x845 mm	730x1204x968 mm	730x1554x968 mm
Weight (Net/Gross)	34/38 kg	86/90 kg	102/110 kg
Input Power	260 W	265 W	300 W
Power Supply	220V, 50Hz		



Biolab Scientific Limited

91 Trafford Crescent, Markham, Ontario, L3R 7J3, Canada

Email: contact@biolabscientific.com | Phone: +1 707 533 1445 | Website: www.biolabscientific.com