



PRODUCT CATALOG



LABORATORY HORIZONTAL AUTOCLAVE



www.biolabscientific.com

300 SERIES HORIZONTAL AUTOCLAVE



High Quality Stainless steel chamber with drying function
 Safety features for over temperature and pressure
 Equipped with a built in printer
 Automatic opening of safety valve if pressure exceeds 0.24MPa
 Door remains locked if pressure exceeds 0.027MPa
 Over temperature, over pressure, auto protection device
 Automatic power cut off and water cut off in case of water shortage

LABORATORY HORIZONTAL AUTOCLAVE

SPECIFICATIONS

Model	BAHZ-301	BAHZ-302	BAHZ-303
Capacity	150 L	200 L	280 L
Adjustment of Temperature	40°C-134°C		
Sterilization Temperature	134°C		
Sterilization Time	0-99 min		
Heat Average	±2°C		
Working Pressure	0.22 MPa		
Time of Drying	0-99 min		
Chamber Dimension	480×850 mm	515×1000 mm	600×1000 mm
Overall Dimension (LxWxH)	1400x600x1300 mm	1400x670x1650 mm	1400x770x1780 mm
Packing Size	800 x 1650 x 1600 mm	860 x 1650 x 1900 mm	960 x 1650 x 2000 mm
Weight	270/360 kg	290/380 kg	340/440 kg
Power	6000 W	9000 W	12000 W
Power Supply	380V, 50Hz	380 V, 50Hz	