



GENERAL PURPOSE INCUBATOR BIGP-204

GENERAL PURPOSE INCUBATOR BIGP-204

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.

Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-204 GENERAL PURPOSE INCUBATOR



Microprocessor controller (with Timer function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely.(Option)

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature.(Option)

SPECIFICATIONS

Model	BIGP-204
Capacity	80 L
Temperature Range	RT +5-65°C
Ambient Temperature	+5-35°C
Temperature Uniformity	±1.5°C (@37°C)
Display Resolution	0.1°C /±0.5°C
Timer Range	0-5999 min
Internal Dimension	500Wx400Dx400H mm
Exterior Dimension	780Wx530Dx560H mm
Shelves	2
Power Consumption	400 W
Electrical Requirement	200 V, 50 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com