



## GENERAL PURPOSE OVEN BOGP-102

## GENERAL PURPOSE OVEN BOGP-102

Designed for continuous drying operations , general purpose ovens provide temperature homogeneity regulated by microcomputer control system. Stainless steel interiors with easy to clean surface while silicon rubber sealing prevents heat loss and gives stability during operation.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

## BOGP-102 GENERAL PURPOSE OVEN



- Micro-computer intelligent control system with LED display
- Forced air convection maintains temperature uniformity
- High quality stainless steel inner chamber with electrostatic spray treated exterior
- Large observation window
- Easy to clean semicircular corners with adjustable shelves
- Silicon rubber sealing ring prevents heat loss
- Over temperature alarm with timer

## SPECIFICATIONS

|                            |                 |
|----------------------------|-----------------|
| Model                      | BOGP-102        |
| Capacity                   | 70 L            |
| Temperature Range          | RT+5-250°C      |
| Temperature Accuracy       | ±1°C            |
| Temperature Uniformity     | ±4.5°C at 250°C |
| Temperature Resolution     | ±1°C            |
| Timer Range                | 0-999 min       |
| Internal Dimension (WxDxH) | 420x450x350 mm  |
| Overall Dimension (WxDxH)  | 630x740x550 mm  |
| Package Size               | 750x820x630 mm  |
| Shelves                    | 2               |
| Weight                     | 44.5/47.5 kg    |
| Power                      | 1200 W          |
| Power Supply               | 220V, 60Hz      |

## OPTIONAL ACCESSORIES

| Accessory Code | Name                                        | Description |
|----------------|---------------------------------------------|-------------|
| BX-59444       | Intelligent program temperature controller  |             |
| BX-59456       | Independent limit of temperature controller |             |





Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)